

## Full-Automatic Filter Integrity tester Integtest V1.1

### Model : IntegtestV1.1

Fully automated integrity test instrument, which performs bubble point, diffusion flow, enhanced bubble point, and water-based test for hydrophobic filters, economic and good performance.

### Main features

- Display of test data and curves in real time, monitoring the whole process of testing.
- Provides automatic printing set function, the user can operate more simple and convenient.
- Water Based test (WH) for Hydrophobic Filter: using purified water as test liquid instead of IPA and Ethanol, thus it can avoid trace quantity of ethanol or IPA contamination when testing filters.
- The instrument can store 500 history record and curve.

### Applied range

- Disc Membrane:  $\Phi 25\text{mm}$ - $\Phi 300\text{mm}$ ;
- Standard Cartridge: 2.5"- 40"
- Capsule & Mini Cartridge
- Air filter test : 2.5"- 40"